

dash

3 ⏻ + ☀️

Scan this QR code to download the full user guide, or type the following link into your web browser: <https://www.lumie.com/support/product-guides/Dash>

Safety

 Read instructions for use.

If the unit is damaged, please contact Lumie before use.

 Keep away from water and damp.

 For indoor use only.

The unit gets warm when the light is on, so make sure that it is not covered.

Use only the mains power adaptor supplied.

Keep the mains power adaptor and cable away from heated surfaces. If either are damaged, they should be replaced with an approved Lumie product.

Only to be used by adults and children over 7 years old. Children using the light, or playing near it, should be supervised.

 Caution – do not attempt to service or modify any part of the unit.

If the unit has been stored outside the recommended operating temperature range (5°C to 35°C), allow it to stand for at least 1 hour before applying power to avoid condensation or overheating.

The light source contained in this luminaire shall only be replaced by Lumie or an authorised Lumie distributor.

Electromagnetic compatibility

Electrostatic discharge may affect the working of this product.

If this occurs, simply switch off for at least 10 seconds and the unit will reset.

Disposal

Do not throw out the mains power adaptor or the unit with normal domestic waste. Take them to an official recycling point to be disposed of responsibly.

WARNING!

Check with your doctor before using this product if any of the following apply:

- You have a major eye problem or have had eye surgery;
- You have been advised to avoid bright light or are taking medication that may cause photosensitivity;
- You have been diagnosed with depression or other mood/sleep disorders. This is especially important if you are also taking medication for this, since light therapy can affect your mood.

Warranty

This unit is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If you need to return the product for servicing, use the original packaging (including internal packing) if possible.

Make sure the mains power adaptor is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first.

Please contact Lumie for further advice.

SN Serial number (see base of unit)

Importer

Use of Dash

How to use Dash for bright light therapy?

Sunlight is a natural mood-booster, influencing body chemistry (e.g. melatonin and serotonin levels). If you're not getting enough exposure to daylight, bright light therapy with Dash is the next best thing, putting you in a better mood, boosting energy, and helping you feel more alert. It can also help relieve the symptoms of winter blues. Position Dash at arm's length (about 50 cm) with the light angled towards your face and set the light intensity to full brightness. Although you don't need to stare at the light source, the light must reach your eyes to have any effect. This means you can use Dash whilst reading, working on your laptop, crafting, etc.

How often should you use Dash for wellbeing support?

Your body's natural rhythms, which influence mood, sleep, and general wellbeing, are governed by how much light (natural and artificial) you're exposed to. Dash mimics the colour and intensity of natural light, so if you're indoors, a bright light therapy session works quickly to improve mood, energy, focus and performance. Used daily, Dash can encourage a more regular sleep/wake pattern, making it easier to get to sleep and wake up feeling more refreshed.

How often should you use Dash for winter blues?

Look out for early signs of winter blues in Sept/Oct and if you start feeling more tired, anxious or irritable, start using Dash. As spring comes around, you can cut back on your bright light sessions. If symptoms return, simply use Dash again for a few days.

When should you use Dash?

If you're using Dash to help with the winter blues, bright light in the morning is most effective, so try that first. For some people, early evening light exposure is beneficial. If you find it difficult to get going in the morning, switch on Dash as soon as you can. If you don't have much time in the morning, then you should aim for a top-up or full session in the afternoon. If you tend to feel sleepy or fall asleep too early, use Dash in the late afternoon/early evening. It's best to avoid bright light too close to bedtime - within about three hours - as you may find it difficult to get to sleep.

How long should you use Dash for?

With the light at full brightness, you're getting the maximum benefit, so position Dash at arm's length (approx. 50 cm) and start with 30 minutes a day (1,500 lux). You may find this is enough, but if you don't notice much of an effect, increase your exposure each day until you find something that works for you. If you're more comfortable with Dash positioned further away, less light will reach your eyes, and you'll need to use it for longer. Similarly, if you're using Dash at reduced brightness, less light will reach your eyes, and you'll need to use it for longer.

There's no harm in using the light for longer than suggested. You can get your bright light therapy all in one go or use Dash on-and-off throughout the day for a cumulative effect. If you find your eyes ache slightly, try using it further away, reduce the brightness, or switch off for a while. Whether you're using Dash for winter blues or for general wellbeing, once you've found a pattern that's right for you, try to use it around the same time every day.

Technical description

Subject to change without notice.

Lumie Dash

Detachable accessory:	Mains power adaptor (supplied) <input type="checkbox"/>
Rated Input:	100–240 V AC, 50/60 Hz, 0.2 A max. (PF > 0.9)
Rated Output:	12.6 V \Rightarrow 0.85 A (12.0 W)

Operating conditions	Temperature:	5°C to 35°C
	Relative humidity:	15% to 93% (non-condensing)
	Atmospheric pressure:	70 kPa to 106 kPa

Storage/transport conditions	Temperature:	-25°C to +70°C
	Relative humidity:	15% to 93% (non-condensing)
	Atmospheric pressure:	25 kPa to 106 kPa
	Colour temperature:	4000 K
	Colour rendering index:	80

At maximum brightness, Dash delivers 10 000 lux at 12 cm.

This product contains a light source of energy efficiency class D.

 European Declaration of Conformity

 UK Conformity Assessed

Contact us

If you have any questions or problems regarding this product, or have suggestions that you feel could help us improve our products or service please email info@lumie.com.

 Outside In (Cambridge) Limited, 3 The Links,
Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:
Outside In (Cambridge) Limited, incorporated in England and Wales.
Registered Number: 2647359.

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© Lumie 2023
LUMIE is a registered trademark of
Outside In (Cambridge) Limited, registrerat i England och Wales.
Outside In (Cambridge) Ltd.
of Outside In (Cambridge) Limited.